ESPAÑOL

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11. Etiquetado

2.6.1 Cutting instruments / indications

Instruments designed to cut different types of tissue. It is recommended to use these instruments carefully and precisely, and ensuring a firm grip on them.

⚠️ A visual inspection of the cutting edge should be made before use. In case of significant deterioration, the instrument should not be used in the intervention, due to the possible loss of precision when making the cut.

2.6.2 Cutting instruments / contraindications

Avoid handling or touching the edge of these instruments to prevent deterioration and possible injury to the user.

⚠️ Do not hit these instruments under any circumstances, they are not designed for this purpose, and since they have a sharp work area, the impact could mean losing control of it, putting the user and the patient at risk.

⚠️ For the same reason explained above, do not lever these instruments under any circumstances, the user or the patient can be injured very easily.

2.7.1 Guidance instruments / indications

Instruments whose main function is to guide other instruments where they are required in the surgical field, or to fix guiding instruments in their correct position.

Many of these instruments have a sharp tip for making pre-drilled holes, or for temporarily penetrating tissues (less than 1 hour of use). It is important to handle them carefully and precisely to ensure correct guidance during the intervention.

Drill guides can be driven into the target tissue by impact, or twisting. In both cases, it is recommended to carry out this process carefully, always checking that the position of the guide is the correct one.

⚠️ A visual inspection of the cutting edge should be made before use. If significant deterioration is detected, the instrument should not be used in the intervention, due to possible loss of precision when drilling.

2.7.2 Guiding instruments / contraindications

Avoid handling or touching the tip of these instruments to prevent premature deterioration. Avoid prying with these instruments, they are not designed for this purpose and can be easily damaged.

If the guide drills by impacting, it is important to impact carefully, precisely and on the same axis as the guide to avoid bending in the process.

2.8.1 Adapters / indications

Instruments that serve to join or adapt an instrument to a surface, or two or more instruments to each other, usually to adapt from one type of anchorage to a different one.

⚠️ It is important to ensure that, once the adaptation has been carried out, the instruments have been firmly secured to each other.

⚠️ It is important to test the adapters before surgery to ensure that they function correctly. If any deterioration is detected, their participation in the intervention must be suspended.

2.8.2 Adapters / contraindications

Avoid hitting the adapters directly. In any case, if the adapted instruments are for impact, the blow must be made on the instruments, and the adapter must only be a transmitter of that force.

Avoid prying the adapter. It is especially important not to pry while the adapter is in use together with your instruments, as this may deteriorate it, and cause damage to other instruments as well.

2.9.1 Trays / indications

Containers to handle and transport other surgical instruments quickly and safely.

It is recommended to handle them carefully to avoid premature deterioration, avoiding blows to both their surface and their interior.

2.9.2 Trays / contraindications

The containers should not be used to store instruments for which they have not been designed.

An inspection of the condition of the container must be carried out to ensure its correct condition.

3. General Warnings

3.1 Scope of use

Valmed SL instruments must be used solely and exclusively in the surgical field. Their use outside of this scope, as well as their misuse and care, can lead to their premature deterioration.

3.2 Allergies or incompatibilities

The instruments must not be used under any circumstances if the user is aware of any type of allergy or incompatibility that the patient may have with the materials with which the Valmed instruments are made. These are listed below:

- AISI 304, AISI 630, AISI 420, AISI 440B, AISI 316LVM stainless steels

- Delrin POM-C and Tecafom technical plastics.

(If you need to know the specific material with which your instrument is made, you can contact the company and it will be communicated to you)

3.3 Surface condition

Valmed instruments will show signs of corrosion if they are exposed to aggressive substances for prolonged periods of time, so it is essential to strictly follow the processing and sterilization instructions.

3.4 Condition of use

For instruments to be in optimal condition, they must be checked and tested prior to use in the operating room.

3.5 Storage

We recommend storing the instruments in dry and clean places.

3.6 Presentation of the instruments

⚠️ The products manufactured and marketed by Valmed are sent to customers in a NON-STERILE state, and the duty to clean and sterilize them before putting them into operation relays on them, following the instructions in section 6. These instructions must be followed whenever the use of the instrument in the surgical field is required.

4. Reparations and responsibilities

Valmed offers a maintenance and repair service for all types of instruments, whether its own instruments or those of third-party companies. However, it reserves the exclusive right to be the only company qualified to carry out the repair and modification of its own instruments. Any repair or modification of its instruments carried out by third parties or companies entails the loss of the product warranty. In addition, Valmed will not be held responsible for possible damages that your instrument may cause from that moment on.

Valmed SL is not responsible for the damage suffered by its instruments when these are caused by improper use or maintenance. This applies particularly in the event of any use that occurs outside the purpose of use stated by the company. It also applies when a repair is carried out that restores or changes the instrument, when it is carried out by companies other than Valmed.

5. Shelf life

The useful life of the products essentially depends on the correct use, care and processing given to them. The guarantee offered by Valmed on its instruments, provided that the aforementioned conditions are met, is 1 year.

After this period, the user is responsible for checking if the instrument is still suitable for use, or if it requires any preventive or corrective maintenance. If necessary, they must always be carried out at the Valmed SL facilities.

6. Processed

Note. Other processes may be suitable, as long as they follow standardized guidelines for the processing of reusable surgical material.

6.1 During use

The correct processing of the instruments begins during the surgical intervention itself. Once the instrument has performed its function, it must be cleaned in the best possible way to prevent traces of aggressive substances such as blood, medications, etc. from remaining dry on its surface. This simple gesture greatly facilitates the task of preserving the instrument.

It is also important to avoid long waiting times between processing steps to prevent these remains from hardening.

6.2 Preparation for decontamination

If possible, instruments should be disassembled prior to further processing steps.

It is important to rinse and clean hard to reach areas such as blind holes, corners and recesses.

The products must be placed in trays of adequate size, and within them, the instruments must keep a minimum distance from each other, it being very important to prevent them from touching or overlapping, to avoid damage and ensure complete cleaning.

6.3 Cleaning and disinfection

Valmed recommends:

A pre-wash with water for 1 minute, followed by a wash with water containing an alkaline cleaning agent at 55°C for 5 minutes, followed by a neutralization with water for 3 minutes, finally rinsing with more water for 2 minutes. After this, drying is carried out at 60° for at least 30 minutes. If necessary, it will be finished by drying with a cloth or with compressed air.

However, we strongly recommend following the instructions provided by the manufacturers of the washing and disinfection machines. In case of conflict with our recommendations, those of the manufacturer of the machines prevail.

6.3 Sterilization

Sterilization must be carried out with devices and packaging materials that comply with current regulations (eg: DIN, EN, ISO).

The autoclave manufacturer's instructions should always be followed.

Valmed recommends for thermally stable products such as instruments saturated steam at 134°, 3 bar, 5 minute plateau time.

The quality of the water and steam that feeds the boiler must be in accordance with what is described in EN standard 285

⚠️ The person in charge of processing the instruments is responsible for ensuring that the process is carried out properly, following globally accepted standards, and that it provides satisfactory results, using the appropriate equipment for it. This usually requires a validation process and routine monitoring of the process to ensure its quality.

7. Functionality test

Instruments should be routinely tested before each use. Special attention should be paid:

- To its general condition (cleanliness, absence of stains, debris, corrosion, discoloration)
- To its operation (ensure that the moving parts behave normally, that the sharp instruments have not lost their sharpness, that the locking elements function normally, that the handles are in good condition, etc.)
- That no component of the instrument is missing (loss of screws, nuts, pins, and other components)
- That they do not present damage (bumps, cracks, bent areas, degraded threads, broken parts, etc.)

We do not recommend using instruments that have lost their marks, due to the loss of traceability that this could entail. Instruments with considerable damage should also not be used. Replacement or sending to technical service for repair should be considered.

8. Storage

Instruments must be stored in dry, clean, dust-free places with a constant temperature. They should be stored in individual packages that offer protection.

During handling and transportation, they must be handled with care to avoid damage.

If they have been stored in a sterile condition, it must be ensured that their sterile condition remains correct when they are to be used.

9. Product Disposal

Valmed does not provide specific instructions for the disposal of its products, since they are not made up of materials that require special handling. The regulations applicable at the national level must be followed.

10. Pictograms used (EN 980)

PICTOGRAM	MEANING	PICTOGRAM	MEANING
	Warning, must be read to avoid harm		Manufacturer symbol, with its information
	Product is sell in a NON STERILE state		Batch number
	Instructions must be read		Batch number, manufacture day
	The product can be sterilized with autoclave at 134°		

11. Labeling