Calle Villa de Bilbao 40,  
46988, Pol. Ind. Fte. Jarro, Paterna, Valencia  
Telf. 96 132 24 20  
gestion@ivalmed.com  
www.ivalmed.com

FR-037 Ed. 5 del 26/02/2024

## ESPAÑOL

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## ENGLISH

### 2.6.1 Cutting devices / indications

Devices designed to cut different types of tissue. It is recommended to use these devices carefully and precisely, and ensuring a firm grip on them.

A visual inspection of the cutting edge should be made before use. In case of significant deterioration, the instrument should not be used in the intervention, due to the possible loss of precision when making the cut.

### 2.6.2 Cutting devices / contraindications

Avoid handling or touching the edge of these devices to prevent deterioration and possible injury to the user.

Do not hit these devices under any circumstances, they are not designed for this purpose, and since they have a sharp work area, the impact could mean losing control of it, putting the user and the patient at risk.

For the same reason explained above, do not lever these devices under any circumstances, the user or the patient can be injured very easily.

### 2.7.1 Guidance devices / indications

Devices whose main function is to guide other devices where they are required in the surgical field, or to fix guiding devices in their correct position.

Many of these devices have a sharp tip for making pre-drilled holes, or for temporarily penetrating tissues (less than 1 hour of use). It is important to handle them carefully and precisely to ensure correct guidance during the intervention.

Drill guides can be driven into the target tissue by impact, or twisting. In both cases, it is recommended to carry out this process carefully, always checking that the position of the guide is the desired one.

A visual inspection of the cutting edge should be made before use. If significant deterioration is detected, the instrument should not be used in the intervention, due to possible loss of precision when drilling.

### 2.7.2 Guiding devices / contraindications

Avoid handling or touching the tip of these devices to prevent premature deterioration.

Avoid prying with these devices, they are not designed for this purpose and can be easily damaged.

If the guide drills by impacting, it is important to impact carefully, precisely and on the same axis as the guide to avoid bending in the process.

### 2.8.1 Adapters / indications

Devices that serve to join or adapt an instrument to a surface, or two or more devices to each other, usually to adapt from one type of anchorage to a different one.

It is important to ensure that, once the adaptation has been carried out, the devices have been firmly secured to each other.

It is important to test the adapters before surgery to ensure that they function correctly. If any deterioration is detected, their participation in the intervention must be suspended.

### 2.8.2 Adapters / contraindications

Avoid hitting the adapters directly. In any case, if the adapted devices are for impact, the blow must be made on the devices, and the adapter must only be a transmitter of that force.

Avoid prying the adapter. It is especially important not to pry while the adapter is in use together with your devices, as this may deteriorate it, and cause damage to other devices as well.

### 2.9.1 Trays / indications

Containers to handle and transport other surgical devices quickly and safely.

It is recommended to handle them carefully to avoid premature deterioration, avoiding blows to both their surface and their interior.

### 2.9.2 Trays / contraindications

The containers should not be used to store devices for which they have not been designed.

An inspection of the condition of the container must be carried out to ensure its correct condition.

## 3. General Warnings

### 3.1 Scope of use

Ivalmed SL devices must be used solely and exclusively in the surgical field. Their use outside of this scope, as well as their misuse and care, can lead to their premature deterioration.

### 3.2 Allergies or incompatibilities

The devices must not be used under any circumstances if the user is aware of any type of allergy or incompatibility that the patient may have with the materials with which the Ivalmed devices are made. These are listed below:

-AISI 304, AISI 630, AISI 420, AISI 440B, AISI 316LVM stainless steels

-Delrin POM-C and Tefalorm technical plastics.

(If you need to know the specific material with which your instrument is made, you can contact the company and it will be communicated to you)

### 3.3 Surface condition

Ivalmed devices will show signs of corrosion if they are exposed to aggressive substances for prolonged periods of time, so it is essential to strictly follow the processing and sterilization instructions.

### 3.4 Condition of use

For devices to be in optimal condition, they must be checked and tested prior to use in the operating room.

### 3.5 Storage

We recommend storing the devices in dry and clean places.

### 3.6 Presentation of the devices

The products manufactured and marketed by Ivalmed are sent to customers in a NON-STERILE state, and the duty to clean and sterilize them before putting them into operation rests on them, following the instructions in section 6. These instructions must be followed whenever the use of the instrument in the surgical field is required.

## 4. Servicing and responsibilities

Ivalmed offers a maintenance and repair service for all types of devices, whether its own devices or those of third-party companies. However, it reserves the exclusive right to be the only company qualified to carry out the repair and modification of its own devices. Any repair or modification of its devices carried out by third parties or companies entails the loss of the product warranty. In addition, Ivalmed will not be held responsible for possible damages that your instrument may cause from that moment on.

Ivalmed SL is not responsible for the damage suffered by its devices when these are caused by improper use or maintenance. This applies particularly in the event of any damage that occurs outside the purpose of use stated by the company. It also applies when a repair is carried out that restores or changes the instrument, when it is carried out by companies other than Ivalmed.

## 5. Lifetime

The useful life of the products essentially depends on the correct use, care and processing given to them. The guarantee offered by Ivalmed on its devices, provided that the aforementioned conditions are met, is 3 years.

After this period, the user is responsible for checking if the instrument is still suitable for use, or if it requires any preventive or corrective maintenance. If necessary, they must always be carried out at the Ivalmed SL facilities.

## 6. Processing

Note. Other processes may be suitable, as long as they follow standardized guidelines for the processing of reusable surgical material.

### 6.1 During use

The correct processing of the devices begins during the surgical intervention itself. Once the instrument has performed its function, it must be cleaned in the best possible way to prevent remains of aggressive substances such as blood, medications, etc., from remaining dry on its surface. This simple gesture greatly facilitates the work of preserving the instrument. Blind holes, threads and complicated geometries must be washed specifically, using a compressed water gun at 3 bars for at least 15 seconds to remove remains.

It is also important to avoid long waiting times between processing steps to avoid hardening of these materials, a maximum of 2 hours is recommended.

### 6.2 Preparation for decontamination

If possible, instruments should be disassembled before the next processing steps.

It is important to rinse and clean hard-to-reach areas, such as blind holes, threads, corners and recesses.

The products must be placed in trays of appropriate size, and within them, the instruments must keep a minimum distance from each other, it being very important to prevent them from touching or overlapping, to avoid damage and ensure complete cleaning.

### 6.3 Cleaning and disinfection

Ivalmed recommends for cleaning your instruments:

Use an automatic washing and disinfection machine whose effectiveness is verified, through its CE marking in accordance with ISO 15883. This must:

- Use sterile water (10 bacteria/ml) with low endotoxins (maximum 0.25 endotoxins unit/ml) for rinsing, such as purified or highly purified water.
- Use filtered air for the drying process (low content of bacteria, particles)

The parameters and cycle recommended by Ivalmed are as follows:

- A pre-wash with water for 1 minute, uncontrolled temperature, <30°
- A wash with water containing an alkaline cleaning agent at 55°C for 10 minutes. For the concentration of the alkaline agent, follow the manufacturer's instructions.
- Neutralization with water for 5 minutes, uncontrolled temperature <30°.
- Thermal disinfection at 95° for 5 minutes.
- Drying at 100° for 30 minutes.

### 6.4 Sterilization

Sterilization must be carried out with devices and packaging materials that comply with current regulations (e.g.: DIN, EN, ISO)

Ivalmed recommends for thermally stable products such as medical devices, saturated steam at 134°, 3 bar, 5 minutes plateau time.

The quality of the water and steam that feeds the boiler must be in accordance with what is described in the EN 285 standard.

The person in charge of processing the devices is responsible for ensuring that the process is carried out properly, following globally accepted standards, and that it provides satisfactory results, using the appropriate equipment for it. This usually requires a validation process and routine monitoring of the process to ensure its quality.

## 7. Functionality test

Devices should be routinely tested before each use. Special attention should be paid:

- To its general condition (cleanliness, absence of stains, debris, corrosion, discoloration)
- To its operation (ensure that the moving parts behave normally, that the sharp devices have not lost their sharpness, that the locking elements function normally, that the handles are in good condition, etc.)
- That no component of the instrument is missing (loss of screws, nuts, pins, and other components)
- That they do not present damage (bumps, cracks, bent areas, degraded threads, broken parts, etc.)

We do not recommend using devices that have lost their marks, due to the loss of traceability that this could entail. Devices with considerable damage should also not be used. Replacement or sending to technical service for repair should be considered.

## 8. Storage

Devices must be stored in dry, clean, dust-free places with a constant temperature. They should be stored in individual packages that offer protection.

During handling and transportation, they must be handled with care to avoid damage.

If they have been stored in a sterile condition, it must be ensured that their sterile condition remains correct when they are to be used.

## 9. Product Disposal

Ivalmed does not provide specific instructions for the disposal of its products, since they are not made up of materials that require special handling. The regulations applicable at the national level must be followed.

## 10. Pictograms used (EN 980)

| PICTOGRAM | MEANING  | PICTOGRAM | MEANING                                   |
|-----------|--|-----------|---|
|           | Warning, must be read to avoid harm                  |           | Manufacturer symbol, with its information |
|           | Product is sold in a NON-STERILE state               |           | Batch number                              |
|           | Instructions must be read                            |           | Batch number manufacture day              |
|           | The product can be sterilized with autoclave at 134° |           |   |

## 11. Labeling

## Products

These instructions for use apply to the products listed in this document.

## Important Notes

Please keep these instructions in a suitable place, that prevents their deterioration, and that they are easily accessible to any user of the surgical devices produced by Ivalmed SL.

Please pay attention to the following symbol. You will find it throughout this document. This indicates important annotations to be taken into account by the user. Failure to follow these directions could result in damage to the instrument, the user, or the patient.

## 1. Scope of Application

The listed products must be used exclusively by duly qualified personnel, and must only be used in the surgical field, in the manner in which it is explained, and in a suitable hygienic environment. It is mandatory that the personnel assigned to use these devices become accustomed to handling them and familiarize themselves with them before use.

## 2. Intended Use

### General

All listed devices are intended for general surgical use. It is important to handle all devices with care, avoiding them falling to the ground, or accidentally hitting other surfaces or devices.

Ivalmed groups its products by family, depending on their use. The indications and contraindications for each of them are found in the following points.

### 2.1.1 Impact devices / indications

Blunt devices, whose main purpose is to impact or extract other devices. They have robust bodies, to withstand the impact forces to which they are exposed.

To make the impact, use the parts enabled for this purpose, such as dollies or hammers. Exercise extreme caution when impacting the dolly, it is important to hit the most centered possible on the metal part of the dolly, to avoid damage to the plastic handles, and to transmit the force in the most appropriate way possible.

### 2.1.2 Impact devices / contraindications

Avoid hitting the sides, or places that can bend, or that are not enabled for that purpose. Avoid prying, its shape maximizes resistance to forces whose direction is the same as the main axis, not lateral forces.

Exercise extreme caution if the instrument has moving parts to carry out the impact, such as sliding hammers. The user must avoid placing his hands in the movement zones to avoid being accidentally hit by the moving part. You should also avoid placing your hands in places where they can be hit, such as dollies.

### 2.2.1 Drilling devices / indications

Devices designed to perforate hard or soft tissues, they have sharp tips adapted for this purpose. The perforation must be carried out with care not to cause damage to the instrument, or to unwanted tissues. It is convenient to ensure the position of the instrument before starting the drilling.

Avoid using these devices if there is significant deterioration in the area being drilled, as this can cause a loss of precision.

### 2.2.2 Drilling devices / contraindications

Avoid hitting the instrument in places not enabled for it. Avoid prying, these devices are not designed for this purpose and can be damaged.

### 2.3.1 Torque devices / indications

Devices designed to facilitate the twisting of other components, such as screws, nuts, or other devices. Conceived to offer a good grip during use, a firm grip must be ensured by the user during the entire twisting process.

If during the twisting process, the user notices that excessive force must be exerted, extreme caution must be exercised to avoid unwanted slippage. It is good practice, if this situation is detected, to stop the torsion and analyze the causes that are causing this excess resistance, before continuing.

### 2.3.2 Torsion devices / contraindications

Avoid hitting the instrument unless it has places enabled for this purpose, such as dollies. Avoid levering in any other direction than torsion, its design maximizes resistance in this last direction, but not in the others, and the instrument may be damaged.

### 2.4.1 Traction devices / indications

Robust devices to move resistant tissues, such as skin and muscles, by levering on a stable surface, and can also pull other devices. They may or may not have moving parts.

Before pulling, it is necessary to ensure that the grip of what is to be pulled is correct and firm, otherwise, sudden loss of grip may occur that may lead to the unwanted displacement of the instrument, and, due to the nature of the applied forces, concussion may be caused to the user or patient.

If there are moving parts, such as zippers with a trigger, the user must check the proper functioning of these elements before surgery. In particular, you must ensure that the instrument is completely stationary in the locked position. If any anomaly is detected, it is imperative not to proceed with its use, due to the risk of sudden loss of grip during surgery.

### 2.4.2 Traction devices / contraindications

Avoid hitting the instrument unless you have places enabled for it, such as dollies.

Do not use the instrument aggressively, it is important, during traction, to gradually apply force until the desired traction is achieved, minimizing the risk of loss of grip. Similarly, it is important to reduce the applied force smoothly and gradually once the traction is to be completed to avoid unwanted movement of the instrument.

### 2.5.1 Grasping devices / indications

Devices whose purpose is to provide a firm grip on different tissues or other devices.

Carry out a visual inspection of the working area of the instrument before use. If a significant deterioration is detected in the work area, it should not be used, since this deterioration could lead to a sudden loss of grip during the intervention.

### 2.5.2 Grasping devices / contraindications

Avoid hitting the instrument, unless it has areas enabled for this purpose, such as dollies. Avoid prying, these devices are not robust enough and could be damaged in the process.